

## Certificate of Analysis

eShop Dummy Customer CH  
used in eShop only  
CH-

**Print Date:** 25-Jan-2020

---

**Product Name:** Tris Buffer, 50mM Soln For  
Use w/ LAL, 30 ml vial  
**Material Number:** S50-642  
**Batch No:** 0000786697  
**Manufacturing Date:** 12-Feb-2019  
**Expiration Date:** 12-Feb-2022

---

<i>Test</i>	<i>RESULT</i>	<i>SPECIFICATION</i>		<i>UNIT</i>
		<i>MIN</i>	<i>MAX</i>	
Sterility Test	Negative			
Endotoxin Content (EU/ml) 1 EU/ml = 1 IU/ml	< 0.005	***	< 0.005	
Functionality	Pass	Pass	***	
pH Test	7.2	7.0	7.4	

**Additional Information:**

### Leah Baltos

Electronically signed by Leah Baltos

Date: 22-MAR-2019 12:33:15 EST

RELEASE ( Inspection Lot: Usage Decision )

Lonza Walkersville Inc.  
8830 Biggs Ford Road  
Walkersville, MD 21793 8415  
Tel (301) 898 7025  
Fax (301) 845 4024