



## Certificate of Analysis Endotoxin Indicator Vial

Product Code: EIV-0001 (Vial) / EIV-025 (Kit)  
Product Description: Endotoxin Indicator Vial  
Lot Number: EV021901  
Production Date: 27 MAR 2019  
Expiration Date: 2021.02

<u>Assay:</u>	<u>Specification:</u>	<u>Result:</u>
Endotoxin Potency (EU per Vial)	>1,000 EU/vial	Pass

### Reference

United States Pharmacopeia (USP), 2012, General Chapter <85>, Bacterial Endotoxins Test (BET)

This is to certify that this product lot was manufactured according to current cGMP Standards, meets FUJIFILM Wako Chemicals U.S.A. Corporation – LAL Division’s product specifications, and complies with all applicable requirements for product release.

Created By: Jason Holliday Date: 07/31/2019

Quality Assurance Approval By: Josh Wald Date: 7/31/2019