

Certificate of Analysis

eShop Dummy Customer CH
used in eShop only
CH-

Print Date: 17-Dec-2019

Product Name: Pyrosperser Kit Dispersing
Reagent
Material Number: N188
Batch No: 0000756917
Manufacturing Date: 22-Oct-2018
Expiration Date: 06-Sep-2020

<i>Test</i>	<i>RESULT</i>	<i>SPECIFICATION</i>		<i>UNIT</i>
		<i>MIN</i>	<i>MAX</i>	
Vial Part Number: F188				
Vial Lot Number	0000742495			
Final Release Test:				
Appearance of Solution	Pass	PASS	***	
Sterility	Negative	Negative	***	
pH Test	8.3	6.0	9.0	
Endotoxin Content (EU/ml)	< 0.010	***	< 0.025	
1 EU/ml = 1 IU/ml				
Functionality	Pass	PASS	***	
Specific Gravity at 20 C	1.2	1.1	1.2	

Additional Information:

A 2% PYROSPERSE solution is used for the endotoxin content and functionality tests.

Raj Patel

Electronically signed by Raj Patel

Date: 22-OCT-2018 16:15:41 EST

RELEASE (Inspection Lot: Usage Decision)

Lonza Walkersville Inc.
8830 Biggs Ford Road
Walkersville, MD 21793 8415
Tel (301) 898 7025
Fax (301) 845 4024