
CERTIFICATE OF ANALYSIS

Product Code: N201
Product: Pyrogen-Free Test Tubes
10 x 75 mm w/Screw Caps

Lot Number: 0000706689
Manufacture Date: 12-Apr-2018
Expiration Date: 31-Mar-2024

TEST (Method)

SPECIFICATIONS

	Min.	Max.	Results
Sterility Test	***	***	Negative
Appearance	PASS	***	Pass
Endotoxin (EU/tube)	***	< 0.005	< 0,005
1 EU/Tube = 1 IU/Tube			
Functionality	PASS	PASS	Pass

This lot has been reviewed by Quality Assurance in compliance with requirements of Lonza's Quality System.

This document was generated from a validated Part 11-compliant electronic system and thus handwritten signatures are not required.